
## ADAPTACIÓN LENTES TÓRICAS

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## TORIC CONTACT LENS FITTING

Too often, people mistakenly believe they can not wear contact lenses because they have astigmatism. At present the high level of Servilens technology makes it possible to manufacture excellent contact lenses for correcting any astigmatism.

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**European Population**  
**32% Astigmatism**  
**35% Myopia**  
**32% Hyperopes**

### SOFT TORIC CONTACT LENS

Servilens soft toric contact lens is stabilized by adding ballast to the front surface, maintaining the cylinder in the base curve, and will correct all types of corneal and residual astigmatism.

Soft toric contact lenses are made from a 55% and 58% 60% water content material, which allows flexibility as well as high permeability. A line is engraved at 6 o'clock for the identification as well as calculating alignment of the lens.

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- Confirm the best axis ( $5^\circ$  or  $10^\circ$ )

3.-Corneal evaluation. Tear film analysis.

4.-Corneal measurements must be taken : For example: R-Horizontal = 7.85- R-Vertical = 7.35 and Axis -  $120^\circ$

5.-Base Curve:  $7.60 + 1.0 = 8.60\text{mm}$

K average + 1.0mm B.C. = 8.60mm

Plano or Plus spheric power, round off down and Minus spheric power, round off up= $8.70\text{mm}$ .

6.-Using the spectacle Rx.  $-4.50 / -3.75 \times 120^\circ$ (Cylinder - always to be minus)

7.-Order the contact lens, correcting for the back vertex(sphere= $-4.50$  and sphere plus cylinder= $-8.25$ ) Lens 55 Toric Rx  $8.70/ 14.5 -4.25 / -3.25 \times 120^\circ$

Note: ONLY Hyperopes, add  $+0.25$  to  $+0.50$  to the sphere. due to latent hyperopia.

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Customer code number in Servilens to evaluation = 000002

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WEB  
[lens55.com](http://lens55.com)

## TORIC CONTACT LENS FITTING

Servilens recommends its toric contact lens calculator for initial contact lens parameter selection.

If acuity vision CL is not so good as AV spectacles, you may apply NIPO test, Over-refraction and Rule LARS.

There are different methods to determine the amount (degrees) of rotation on a soft toric contact lens when stabilization is not as indicated by the manufacturer (either at 6 o'clock, or at 3-9 o'clock, or the degrees of the refractive axis, etc).

Methods:

- 1- A reticule in our slit lamps to measure the amount of rotation.
- 2- Calculating the amount of rotation as a clock
- 3- Over refraction calculator.
- 4- NIPO method. (pag 96)

Servilens recommends its calculator for determining the amount of rotation on a soft toric contact lens.

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## TORIC CONTACT LENS FITTING

---

IFU Instructions For User- Contact Lens Care

Servilens recommends to deliver and sign the IFU. Once the right contact lens has been chosen, instructions for lens insertion and removal, hygiene, and wearing time will be provided for the user.

### Final Contact Lens Prescription Determination

After the initial contact lens exam, the user requires one or more follow-up visits during the adaptation period to assess lens performance and ensure that the contact lenses are providing optimal comfort and clarity while maintaining healthy eyes. This follow-up is crucial to complete the fitting process and the necessary information required for the final contact lens prescription.

### Ongoing Contact Lens Care

Many contact lens complications do not cause symptoms in the early stages, but professional examination can detect subtle changes in vision and eye health. Early detection of problems can prompt corrective measures before more serious complications may develop.